



MA Moisture Analyzer



1. Halogen test lamp
2. 5 inch touch panel
3. Intuitive user interface make users can quickly familiar with the operational process of equipment, one key to start drying method.
4. 50-180°C test temperature
5. High precision sensor, supply good quality for weighing

| | |
|--------------------|--------------------------------------|
| Model No. | MA110-1 |
| Capacity | 100g |
| Readability(g) | 0.001g |
| Repeatability | 0.02-0.1% ($\geq 3g$) |
| Operating Temp. | 5°C-40°C |
| Pan Size | \varnothing 100mm |
| Heat Source | Halogen Lamp |
| Temp. Range | 50°C - 180°C |
| Moisture Range | 0.00%-100.00% |
| Time Setting | 0-99 minutes, interval 1 minute |
| Heater | 110-220V $\pm 10\%$, 50Hz $\pm 1Hz$ |
| Display | 5 inch touch panel |
| Calibration | External Calibration |
| Interface | RS232 |
| Housing Dimensions | 330x180x140mm |
| Packing Size | 470*310*330mm |
| N.G. Weight | 3kg/4.5kg |