



## Specification

Model	SH-WB-10GDH
Temp Range	Ambient +10°C to 300°C (Ambient +18°F to 572°F)
Temp Controller	Digital P I D Controller
Demension Internal (WxDxH)	250×300×170(mm) 9.84×11.81×6.69(inch)
Demension External (WxDxH)	635×385×345(mm) 25.00×15.15×13.58(inch)
Material Internal	Stainless Steel
Material External	Steel Plate With Powder Heating Coated
Capacity	13 L
Heater Capacity	2 KW
Electrical Requirements 230V ,50/60Hz, 1Φ	9.0 A
Electrical Requirements 120V ,50/60Hz, 1Φ	17.3 A
Weight	21 Kg

## Feature

- Automatic reset function incoming power cut
- Buzzer on when set time arrived
- Electric leakage circuit breaker
- Over temperature protection

## Option

- Heater cover
- Bath Lid



## Specification

Model	SH-WB-20GDH
Temp Range	Ambient +10°C to 300°C (Ambient +18°F to 572°F)
Temp Controller	Digital P I D Controller
Demension Internal (WxDxH)	350×400×170(mm) 13.78x15.75x6.69(inch)
Demension External (WxDxH)	732x485x345(mm) 28.81x19.09x13.58(inch)
Material Internal	Stainless Steel
Material External	Steel Plate With Powder Heating Coated
Capacity	24 L
Heater Capacity	3.5 KW
Electrical Requirements 230V ,50/60Hz, 1Φ	16.0 A
Electrical Requirements 120V ,50/60Hz, 1Φ	30.7 A
Weight	29 Kg

## Feature

- Automatic reset function incoming power cut
- Buzzer on when set time arrived
- Electric leakage circuit breaker
- Over temperature protection

## Option

- Heater cover
- Bath Lid



## Specification

Model	SH-WB-40GDH
Temp Range	Ambient +10°C to 300°C (Ambient +18°F to 572°F)
Temp Controller	Digital P I D Controller
Demension Internal (WxDxH)	600×400×200(mm) 23.62x15.75x7.87(inch)
Demension External (WxDxH)	1035×500×390(mm) 40.74x19.68x15.35(inch)
Material Internal	4.8 KW
Material External	Stainless Steel
Capacity	48 L
Heater Capacity	Steel Plate With Powder Heating Coated
Electrical Requirements 230V,50/60Hz, 1Φ	23.0 A
Electrical Requirements 120V ,50/60Hz, 1Φ	44.1
Weight	38 Kg

## Feature

- Automatic reset function incoming power cut
- Buzzer on when set time arrived
- Electric leakage circuit breaker
- Over temperature protection

## Option

- Heater cover
- Bath Lid