

RH-24 3D Gyrotory Rocker

RH-24 3D Gyrotory Rocker is an economical 3D gyrotory rocker with perfect 3-dimensional motion for gentle but thorough mixing. It is an ideal choice for small sample mixing. It suits for blood collection tube and centrifuge tube, applies for all blot related applications such as gel staining and destaining. The concave pad can firmly keep the tubes. RH-24 is suitable for 1.5mlx60, 3mlx30, 7mlx21, 5ml x5. suitable for cold closet.



【 Product Features 】

1. Gentle but thorough mixing
2. Suits for blood collection tubes and centrifuge tubes
3. Tilt angle and speed have been optimally set
4. Safe for use in cold rooms and incubators

【 Technical Parameters 】

Model	RH-24
Speed Range	20rpm (fixed)
Shaking Orbit	3D / 20° (fixed)
Max. Load	0.8kgs
Voltage	AC 220V / AC 110V, 50/60Hz
Power	45W
Dimension (WxDxH)	240 x 170 x 150mm
Net Weight	1.6kgs



HS-25/HS-25A Rocking Shaker

HS-25 Rocking Shaker is widely used in immobilization of running gel, shaking for dyeing or fading the Coomassie blue, immobilization, dyeing and developing of silver-staining, developing and fixation of the X-ray film in autoradiography experiment, and further processing of cellulose membrane after gel shift such as molecular hybridization, staining with antigen-antibody. It is also used for cell culture and cell membrane transfer.

【 Product Features 】

1. Platform made of stainless steel and natural rubber, which is alkali-resistant, non-deformation, easy-clean and spill-prevention.
2. Speed real-time display and diminishing time display.
3. The bearing is articulated rotating, secure and no wear, suitable for long time use.
4. Compact designed, take very little space.
5. Brushless DC motor, free maintenance.

【 Technical Parameters 】

Model	HS-25	HS-25A
Speed Range	10 ~80rpm	10 ~80rpm
Swing Cycloid	25mm (Up and Down) / 10°	25mm (Up and Down) / 10°
Time Range	1min ~99h59min or continuous	1min ~99h59min or continuous
Max. Load	5.0kgs (single-deck platform)	5.0kgs (double-deck platform)
Platform Dimension	228 x 280mm	228 x 280mm
Voltage	AC 100 ~230V, 50/60 Hz	AC 100 ~230V, 50/60 Hz
Power	50W	50W
Fuse	250V, 1A, Φ5x20	250V, 1A, Φ5x20
Dimension (WxDxH)	284 x 264 x 134mm	284 x 264 x 248mm
Net Weight	2.5kgs	5.2kgs



SciLution

บริษัท ไซลิวชัน จำกัด Scilution Co., Ltd.

Tel : 02 964 8018 | Line ID : @scilution

www.scilution.co.th / www.scilution.net / www.scilution.com