

# Vertical Shaker

## OLT-06A



Vertical Shaker is a kind of laboratory liquid liquid separating & extraction instrument by vertical shaking. Widely used in oil extraction from water, volatile phenol extraction, food & grease extraction, pesticide residue extraction, water organic pollutant extraction, soil harmful substance extraction, etc.

### Specification

Model	OLT-06A
Shake Speed	20-100 HZ
Shake Range	40 mm
Shake Method	Vertical Shake (tilt angle 0-20° adjustable)
Frequency Control	Infinitely variable speed, LCD display
Time Set Range	0-9999 Sec (count down applied)
Applied Vessel	50-1000 mL
Vessel Number	6 pcs
Extraction Efficiency	95-110%
Load Rated	Bilateral symmetry 11kg
Power Supply	220V, 50/60Hz; 350W

