

Washer Disinfector

BK-LW Series



The lab glassware washers integrate programmed cleaning, disinfecting and drying functions to clean all glassware stuff, for beaker, flask, pipette, petri dish, burette, measuring cylinder, test tube etc.

Features:

- ✓ High efficient cleaning system, designed with European reliable pump, optimized spray arm and nozzles, and self-cleaning program;
- ✓ Rapid and efficient drying system, designed with independent air heater, HEPA filter and drying process;
- ✓ Water heating temperature can reach 99 °C, while hot air drying temperature can reach 120 °C;
- ✓ Professional configuration design—stainless steel for corrosion resistance; reinforced glass window mounted to achieve a clear view on washing process;
- ✓ PLC touch screen control with preset 6 programs and 99 customized programs
- ✓ Safety protection—electronic security door lock to prevent unexpected door opening, water & drying air temperature dual control and equipped with emergency switch.

Features:

Model	BK-LW120	BK-LW220	BK-LW320
Capacity	120L	220L	320L
Water Consumption/Cycle		15L	
Inlet Water Pressure		0.03~0.8Mpa	
Noise	≤50dB	≤55dB	
Material	Internal material: SS316, Outer material: SS304		
Overall Power Consumption	6KW	6.5/16.5KW	6.5/18KW
Power Consumption-Water Washing	0.75KW		1.5 KW
Power Consumption-Water Heating	5KW	5/15KW	
Power Consumption-Drying	2.1KW	3.6KW	
Power Supply	AC220V, 50Hz	AC220V, 50Hz; 380V, 50Hz	
Cleaning Racks	1 pcs	2 pcs	3 pcs
Standard Cleaning Racks	BKIR36*1	BKIR36*2	BKIR36*3
Standard Accessories	One bucket RBS A 155 leaner (1 L), One bucket RBS A 375 neutralizer (1 L)		
Optional Accessories	USB interface, printer, cleaning racks		
Washer Chamber Size(W*D*H)mm	600*629*465	600*629*658	600*629*820
External Size(W*D*H)mm	960*740*850	780*690*1760	790*690*1950
Package Size(W*D*H)mm	1100*950*1050	1000*850*2100	1000*840*2180
Gross Weight(kg)	145	330	350

