

MAGNETIC STIRRER (MULTI POSITIONS)

ST-300 series



Stainless Steel #304



Speed Adjustable



Strong Power



Rechargeable Battery

Model	Position	Dimension (cm)	Speed	Capacity
ST-300-1C With battery	3	14 x 55 x 5.5	100 - 1200 rpm	1000 ml
ST-300-3C	3	14 x 55 x 5.5	100 - 1200 rpm	1000 ml
ST-300-3C With Battery	3	14 x 55 x 4.5	100 - 1200 rpm	1000 ml

*Note : Custom-size is Available. (1C = 1 Controller, 3C = 3 Controller)

*Note : Battery life time 4-5 hours

