

# MAGNETIC STIRRER (MULTI POSITIONS)

## ST- series



ST-500



ST-900



ST-1200



Stainless Steel #304



Speed Adjustable



Strong Power



Separate Control

Model	Position	Dimension (cm)	Speed	Capacity
ST-500	5	33 x 36 x 9	100 - 1200 rpm	1000 ml
ST-900	9	40 x 42 x 9	100 - 1200 rpm	1000 ml
ST-1200	12	44 x 35 x 9	100 - 1200 rpm	1000 ml
ST-1500	15	55 x 35 x 9	100 - 1200 rpm	1000 ml
Options	LED : LED Light for any positions			

\*Note : Custom-size is Available.

