

MAGNETIC STIRRER (SINGLE POSITIONS)

ST- series

ST-ECO



ST-GO



Stainless Steel #304



Strong Power



Speed Adjustable



Analog controller

Model	Position	Dimension (cm)	Speed	Capacity
ST-ECO	1	20.5 x 20.5 x 5	100 - 1200 rpm	1000 ml
ST-GO With battery	1	20.5 x 20.5 x 5	100 - 1200 rpm	1000 ml

*Note : Battery life time 4-5 hours



SciLution

บริษัท ไซลิวชัน จำกัด Scilution Co., Ltd.

www.scilution.co.th | Youtube : SCILUTION

Tel : 02 964 8018 | Line ID : @scilution

