

SOUND ROOF BOX

ตู้กันเสียง



SOUND PROOF BOX



 Aluminium Profiler Base

 Acoustic Foam

 Noise Reducer

 Custom Made

Model	Dimension (W x D x H)	Noise Reducer (dB)
For Ultrasonic bath	54 x 53 x 44 cm	50-60
For Sieve Shaker	64 x 84 x 56 cm	50-60
For Tap Density	44 x 64 x 41 cm	50-60

* Custom-size is Available *

*Note : Not 100% Sound Proof

