

# SUS304 Forced Air Drying Oven

## Characteristics:

WGLL-BE series have the following characteristics other than the above functions.

- LCD temperature controller.
- Over-temperature protection system.
- High precision digital independent temperature limiter.



## Main Technical Parameters:

Model	WGLL-30BE (SUS304)	WGLL-45BE (SUS304)	WGLL-65BE (SUS304)	WGLL-85BE (SUS304)	WGLL-125BE (SUS304)	WGLL-230BE (SUS304)	
Cycle Mode	Forced Convection						
Function	Temp. Range						
	Temp. Resolution						
	Temp. Motion						
	Temp. Uniformity						
Structure	Inner Chamber						
	Outer Shell						
	Insulation Layer						
	Heater						
	Power Rating	0.8kW	1.2kW	1.6kW	2.0kW	2.3kW	3.0kW
Exhaust Hole	φ28mm top(with function of test hole)					Rear Exhaust	
Controller	Temp. Control Mode						
	Temp. Setting Mode						
	Temp. Display Mode						
	Timer						
	Operation Function						
	Additional Function						
Sensor	PT100						
Safety Device	High precision digital independent temperature limiter, over temperature sound and light alarm						
Specification	Inner Chamber(W*L*Hmm)	310*310*310	350*350*350	400*360*450	450*420*450	500*450*550	600*500*750
	Exterior Size(W*L*Hmm)	450*500*690	490*540*730	540*550*830	590*610*830	640*640*930	730*680*1250
	Packing Size(W*L*Hmm)	550*585*845	590*625*885	640*635*985	690*695*985	740*725*1085	860*790*1400
	Capacity	30L	45L	65L	85L	125L	230L
	Shelf Number	6	7	9	9	13	17
	Load Per Rack	15kg					
	Shelf Space	35mm					
	Input(50/60Hz)	AC220V/3.6A	AC220V/5.5A	AC220V/7.2A	AC220V/9.0A	AC220V/10.5A	AC220V/13.6A
	NW/GW(kg)	33/37	37/43	44/49	50/56	60/66	94/120
	Accessory	Shelf	2				
Shelf Frame		4					
Optional Accessories	5.0-inch true color touch screen ;USB Data storage						

